



GeneXPlus™ Transfection of Plasmid DNA into TeloHAEC Cells

TeloHAEC (ATCC® Cat. No. CRL-4052™) is an endothelial, normal human aortic cell line. ATCC achieved transfection efficiencies of 45%, using the protocol described below.

General Considerations for using the GeneXPlus transfection reagent:

- All steps should be performed in a biosafety cabinet using proper aseptic technique.
- **Cell conditions.** Cells should be passaged at least once after thaw and the use of low-passage cells is recommended. Passage the cells 18-24 hours before transfection to ensure the cells are actively dividing and that they will be at the appropriate cell density at the time of transfection. Make sure that the cells are healthy and are $\geq 90\%$ viable, prior to transfection.
- **Seeding density.** Cell density should be 50-80% confluent on the day of transfection. See specified seeding density in the individual protocols and in Table 1. *Note: Determine the optimal cell density for each cell type in order to maximize transfection efficiency.*
- **DNA purity.** Use highly purified plasmid preps that are free from phenol or other contaminants. Plasmid DNA preps that are endotoxin-free are desirable.
- **Presence of antibiotics and other inhibitors.** Antibiotics will inhibit transfection complex formation and therefore should be excluded from the complex formation step. Transfection complexes can be added to cells grown in complete culture medium containing serum and low levels of antibiotics if required. Heparin in the medium should be removed at the day of performing the transfection.
- **Complex formation conditions.** Prepare GeneXPlus Reagent and DNA complexes in serum-free growth medium. ATCC recommends using Opti-MEM I Reduced-Serum Medium to dilute the DNA before complex formation.

Materials required:

Material Required	Catalog No.
TeloHAEC cells	ATCC® CRL-4052™
VCBM	ATCC® PCS-100-030™
VEGF KIT	ATCC® PCS-100-041™
FBS	ATCC® 30-2020™
GeneXPlus™	ATCC® ACS-4004
Opti-MEM® I Reduced-Serum Media	Life Technologies™ 31985-062
Plasmid DNA of interest (1µg/µL)	
Tissue culture plates and supplies	

Protocol:



GeneXPlus™ Transfection of Plasmid DNA into TeloHAEC Cells

The following protocol describes how to transfect plasmid DNA into TeloHAEC cells using the GeneXPlus Reagent in a **single well of a 6-well plate**. The reaction may be scaled up as needed. Please refer to Table 1 for recommended reaction conditions for other dish or plate sizes.

A. Preparation of the cells for transfection

The day before transfection:

1. Count and measure cells for density and viability.
2. Plate **5.0×10^4 - 8.0×10^4** cells per well in 2 mL of complete growth medium (ATCC VCBM + VEGF KIT). Cell density should be **50-80%** confluent on the day of transfection.
3. Incubate cells overnight at **37°C** with **5% CO₂**.

The day of transfection:

1. Remove old media.
2. Replace old media with fresh complete growth media **without heparin** to a total volume of 2.0 mL.

B. Preparation of the DNA: GeneXPlus transfection complexes

1. Warm GeneXPlus, plasmid DNA, and Opti-MEM I Reduced-Serum Medium to room temperature and vortex gently to mix.
2. Pipette 200 μ L Opti-MEM I Reduced-Serum Medium into a sterile microcentrifuge tube.
3. Add 2.0 μ L (1.0 μ g/ μ L) plasmid DNA.
4. Mix thoroughly with gently pipetting.
5. Add 4.0 μ L GeneXPlus Reagent to the diluted DNA mixture. *Note: Do not let the pipette tip or the reagent come into contact with the sides of the plastic tube.*
6. Mix GeneXPlus:DNA complexes thoroughly using either a vortex or by pipetting briefly.
7. Collect contents at bottom of the tube using a mini-centrifuge.
8. Incubate GeneXPlus:DNA complexes at room temperature for 15 minutes.

C. Addition of DNA:GeneXPlus transfection complexes to the cells

1. Distribute the complexes to the cells by adding the complexes drop-wise to different areas of the wells.
2. Gently rock the culture vessel back and forth and from side to side to evenly distribute the GeneXPlus:DNA complexes.


D. Post-Transfection Handling

1. Incubate for **24-72** hours. Replace growth media 24 hours post transfection.
2. Wait for 18-24 hours post-transfection before assaying for transgene expression

Table 1: Recommended Reaction Conditions for different size culture vessels.

Culture Vessel	24 well plate	12 well plate	6 well plate	10 cm dish
Surface area	1.9 cm ²	3.8 cm ²	9.6 cm ²	59 cm ²
Complete Growth Medium	0.5 mL	1.0 mL	2.0 mL	10.0 mL
Opti-MEM I Reduced Serum Medium	50 µL	100 µL	200 µL	1.0 mL
DNA (1 µg/µL stock)	0.5 µg	1.0 µg	2.0 µg	10 µg
GeneXPlus Reagent	1.0 µL	2.0 µL	4.0 µL	20 µL

Note: Always include a control condition consisting of an empty vector plasmid or a plasmid expressing GFP.



Transfection efficiency of GeneXPlus reagent on TeloHAEC cells in a 6-well plate. Cells were transfected with pcDNA 6.2 -EF1α-emGFP vector at 2.0 µg DNA with 4.0 µL of reagent (1:2) in Opti-MEM I Reduced Serum Media. No significant signs of cellular toxicity were observed.