

**ATCC medium: 1799 LB medium (ATCC medium 1065) plus 10 mcg/ml neomycin**

ATCC Medium 1065 (see below) plus 10.0 mcg/ml neomycin  
Add sufficient volume of filter-sterilized neomycin to achieve a final concentration of 10 mcg/ml.

*ATCC Medium 1065:*

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.