



THE ESSENTIALS OF
LIFE SCIENCE RESEARCH
GLOBALLY DELIVERED™

Pharmaceutical Microbiology




Table of Contents

Pharmaceutical Microbiology	1
Strains specified by official microbial assays	2
United States Pharmacopeia (USP)	2
European Pharmacopeia (EP) Edition 8.1	5
Japanese Pharmacopeia (JP)	7
Strains listed by genus and species	10

**THE ESSENTIALS
OF LIFE SCIENCE
RESEARCH
GLOBALLY
DELIVERED™**

ATCC provides research and development tools and reagents as well as related biological material management services, consistent with its mission: to acquire, authenticate, preserve, develop, and distribute standard reference microorganisms, cell lines, and related materials for research in the life sciences.

For over 85 years, ATCC has been a leading provider of high-quality biological materials and standards to the life science community. We are an independent, 501(c)(3) non-profit entity focused on scientific enablement at universities, research institutes, government agencies, and commercial research labs. Our diverse and comprehensive resources in cell biology and microbiology have been central to the growth of the biotechnology age. ATCC has as its core mission to source,

authenticate and further develop products and services essential to the needs of basic and applied life science work.

ATCC distributes to more than 165 countries on 6 continents and has a growing international network of 12 distribution partners. Our infrastructure and experience in biological materials logistics enables us to work effectively with researchers no matter where they are located.

PHARMACEUTICAL MICROBIOLOGY

Pharmaceutical companies are well versed in the important role of microbiological testing in its basic functions – product research and development, process validation, manufacturing, and quality control.

The pharmaceutical industry looks to ATCC to provide top-quality, fully characterized strains necessary to maintain the highest levels of product integrity and reputation. With ATCC Genuine Cultures[®], pharmaceutical customers can rely on:

- Polyphasic identification and strain characterization
- Proprietary systems for managing strain stocks
- Meticulous propagation, preservation, and storage protocols
- Robust quality systems
- Over 85 years of service


Whenever possible, ATCC balances traditional biochemical testing methods with automated phenotypic and genotypic analyses. This provides highly accurate identification across a wide array of microorganisms and can help avoid the pitfalls of misidentification, painful recalls, and regulatory repercussions. Each vial contains a minimal passage descendant of the original material and has been handled only by ATCC.

ATCC is a trusted scientific resource and an ISO 9001 certified, ISO 17025 and ISO Guide 34 accredited organization. Our customer and technical service teams, as well as our global network of authorized distributors, understand the complexities of distributing and supporting microbial cultures worldwide. Our laboratory procedures maintain strains carefully, safely, and effectively to ensure optimal viability, backed by our warranty.

Strains provided for key pharmaceutical applications are listed in two ways:

- Strains specified by United States Pharmacopeia (USP), European Pharmacopeia (EP), or Japanese Pharmacopeia (JP) official microbial assays
- Strains listed alphabetically by genus and species paired with relevant applications

To learn more visit our website at www.atcc.org/pharmaceutical



STRAINS SPECIFIED BY OFFICIAL MICROBIAL ASSAYS

Official Microbial Assays: United States Pharmacopeia (USP)

Application	ATCC® No.
USP CHAPTER: <51> ANTIMICROBIAL EFFECTIVENESS TESTING, USP 37-NF32, 2014	
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLR1 034(120)	16404™
<i>Escherichia coli</i> Crooks	8739™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
USP CHAPTER: <61> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS: MICROBIAL ENUMERATION TESTS, USP 37-NF32, 2014	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLR1 034(120)	16404™
USP CHAPTER: <62> MICROBIOLOGICAL EXAMINATION OF NONSTERILE PRODUCTS: TESTS FOR SPECIFIED MICROORGANISMS, USP 37-NF32, 2014	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Escherichia coli</i> Crooks	8739™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium CDC 6516-60	14028™
<i>Candida albicans</i> 3147	10231™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™
<i>Clostridium sporogenes</i>	19404™
USP CHAPTER: <63> MYCOPLASMA TESTS, USP 37-NF32, 2014	
<i>Acholeplasma laidlawii</i> PG8	23206™
<i>Mycoplasma gallisepticum</i> PG 31	19610™
<i>Mycoplasma fermentans</i> PG 18	19989™
<i>Mycoplasma hyorhinis</i> BTS-7	17981™
<i>Mycoplasma orale</i> CH 19299	23714™
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	15531™
<i>Mycoplasma synoviae</i> WVU 1853	25204™
USP CHAPTER: <71> STERILITY TESTS, USP 37-NF32, 2014	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Clostridium sporogenes</i>	19404™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLR1 034(120)	16404™
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
<i>Bacteriodes vulgatus</i>	8482™
USP CHAPTER: <81> ANTIBIOTICS-MICROBIAL ASSAYS, USP 37-NF32, 2014	
<i>Saccharomyces cerevisiae</i>	9763™
<i>Micrococcus luteus</i> 130.21	10240™
<i>Mycobacterium smegmatis</i>	607™
<i>Pseudomonas aeruginosa</i>	25619™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA	29737™
<i>Bordetella bronchiseptica</i> NRRL B-140	4617™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™

Official Microbial Assays: United States Pharmacopeia (USP)

Application	ATCC® No.
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
<i>Saccharomyces cerevisiae</i>	2601™
<i>Enterococcus hirae</i> FDA M19	10541™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> 3R7089 strain Oxford	9144™
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	10031™
<i>Escherichia coli</i> MacLeod	10536™
USP CHAPTER: <91> CALCIUM PANTOTHENATE ASSAY, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
USP CHAPTER: <115> DEXPANTHENOL ASSAY, USP 37-NF32, 2014	
<i>Pediococcus acidilactici</i>	8042™
USP CHAPTER: <171> VITAMIN B12 ACTIVITY ASSAY, USP 37-NF32, 2014	
<i>Lactobacillus leichmannii</i> 313	7830™
USP CHAPTER: <441> NIACIN OR NIACINAMIDE ASSAY, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
USP CHAPTER: <1035> BIOLOGICAL INDICATORS, USP 37-NF32, 2014	
<i>Geobacillus stearothermophilus</i> NCA 26	12980™
USP CHAPTER: <1046> CELLULAR AND TISSUE-BASED PRODUCTS, USP 37-NF32, 2014	
Hybrid Moloney/Amphotropic Murine leukemia virus (Mo/A-MuLv)	VR-1448™
Hybrid Moloney/Amphotropic Murine leukemia virus (Mo/A-MuLv)	VR-1450™
USP CHAPTER: <1072> DISINFECTANTS AND ANTISEPTICS, USP 37-NF32, 2014	
<i>Escherichia coli</i> AMC 198	11229™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Pseudomonas aeruginosa</i> PRD-10	15442™
<i>Candida albicans</i> 3147	10231™
<i>Candida albicans</i> 132	2091™
<i>Penicillium chrysogenum</i> Wis. 49-133	11709™
<i>Aspergillus brasiliensis</i> WLR1 034(120)	16404™
<i>Bacillus subtilis</i> PRD 66	19659™
USP CHAPTER: <1211> STERILIZATION AND STERILITY ASSURANCE OF COMPENDIAL ARTICLES, USP 37-NF32, 2014	
<i>Geobacillus stearothermophilus</i>	7953™
<i>Bacillus atrophaeus</i> NRS 1221A	9372™
USP CHAPTER: <1229.1> STEAM STERILIZATION BY DIRECT CONTACT, USP 37-NF32, 2014	
<i>Geobacillus stearothermophilus</i> NCA 26	12980™
<i>Geobacillus stearothermophilus</i>	7953™
USP DIETARY SUPPLEMENT CHAPTER: <2021> MICROBIAL ENUMERATION TESTS - NUTRITIONAL AND DIETARY SUPPLEMENTS, USP 37-NF32, 2014	
<i>Candida albicans</i> 3147	10231™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium NCTC 74	13311™
<i>Aspergillus brasiliensis</i> WLR1 034(120)	16404™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Escherichia coli</i> Crooks	8739™
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS CAPSULES, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS ORAL SOLUTION, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™

Official Microbial Assays: United States Pharmacopeia (USP)

Application	ATCC® No.
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS TABLETS, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS WITH MINERALS CAPSULES, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS WITH MINERALS ORAL SOLUTION, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: OIL- AND WATER-SOLUBLE VITAMINS WITH MINERALS TABLETS, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
USP DIETARY SUPPLEMENTS: WATER-SOLUBLE VITAMINS CAPSULES, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: WATER-SOLUBLE VITAMINS TABLETS, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
USP DIETARY SUPPLEMENTS: WATER-SOLUBLE VITAMINS WITH MINERALS CAPSULES, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: WATER-SOLUBLE VITAMINS WITH MINERALS ORAL SOLUTION, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
<i>Pediococcus acidilactici</i>	8042™
USP DIETARY SUPPLEMENTS: WATER-SOLUBLE VITAMINS WITH MINERALS TABLETS, USP 37-NF32, 2014	
<i>Lactobacillus plantarum</i> 17-5	8014™
<i>Lactobacillus leichmannii</i> 313	7830™
USP MONOGRAPHS: BIOLOGICAL INDICATOR FOR DRY-HEAT STERILIZATION, PAPER CARRIER, USP 37-NF32, 2014	
<i>Bacillus atrophaeus</i> NRS 1221A	9372™
USP MONOGRAPHS: BIOLOGICAL INDICATOR FOR ETHYLENE OXIDE STERILIZATION, PAPER CARRIER, USP 37-NF32, 2014	
<i>Bacillus atrophaeus</i> NRS 1221A	9372™
USP MONOGRAPHS: BIOLOGICAL INDICATOR FOR STEAM STERILIZATION, PAPER CARRIER, USP 37-NF32, 2014	
<i>Geobacillus stearothermophilus</i> NCA 26	12980™
<i>Geobacillus stearothermophilus</i>	7953™
USP MONOGRAPHS: PENICILLAMINE, USP 37-NF32, 2014	
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
USP MONOGRAPHS: PENICILLIN G PROCAINE AND NOVIOBIOCIN SODIUM INTRAMAMMARY INFUSION, USP 37-NF32, 2014	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> PCI 1209/N	12692™
USP MONOGRAPHS: TETRACYCLINE HYDROCHLORIDE AND NOVIOBIOCIN SODIUM TABLETS, USP 37-NF32, 2014	
<i>Escherichia coli</i> MacLeod	10536™

Official Microbial Assays: European Pharmacopeia (EP) Edition 8.1

Application	ATCC® No.
EP 2.6.1 STERILITY	
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Clostridium sporogenes</i>	19404™
<i>Candida albicans</i> 3147	10231™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Rosenbach FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
EP 2.6.7 MYCOPLASMAS	
<i>Mycoplasma fermentans</i> PG18	19989™
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	15531™
<i>Mycoplasma hyorhinis</i> BTS-7	17981™
<i>Mycoplasma gallisepticum</i> PG31	19610™
<i>Mycoplasma orale</i> CH 19299	23714™
<i>Mycoplasma synoviae</i> WVU 1853	25204™
<i>Acholeplasma laidlawii</i> PG8	23206™
EP 2.6.12 MICROBIAL EXAMINATION OF NON-STERILE PRODUCTS: MICROBIAL ENUMERATION TESTS	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Candida albicans</i> 3147	10231™
EP 2.6.13 MICROBIAL EXAMINATION OF NON-STERILE PRODUCTS: TEST FOR SPECIFIED MICRO-ORGANISMS	
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium CDC 6516-60	14028™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™
<i>Clostridium sporogenes</i>	19404™
<i>Escherichia coli</i> Crooks	8739™
<i>Candida albicans</i> 3147	10231™
EP 2.6.27 MICROBIOLOGICAL CONTROL OF CELLULAR PRODUCTS	
<i>Candida albicans</i> 3147	10231™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™
<i>Propionibacterium acnes</i> Gerath	11827™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Clostridium sporogenes</i>	19404™
<i>Streptococcus pyogenes</i> Bruno	19615™
<i>Bacteroides fragilis</i> VPI 2553	25285™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Rosenbach FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> 33114	9610™
EP 2.6.31 MICROBIOLOGICAL EXAMINATION OF HERBAL MEDICINAL PRODUCTS FOR ORAL USE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> Rosenbach FDA 209	6538™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Escherichia coli</i> Crooks	8739™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium CDC 6516-60	14028™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™

Official Microbial Assays: European Pharmacopeia (EP) Edition 8.1

Application	ATCC® No.
EP 2.7.2 MICROBIOLOGICAL ASSAYS OF ANTIBIOTICS	
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i>	10031™
<i>Micrococcus luteus</i> 130.21	10240™
<i>Escherichia coli</i> MacLeod	10536™
<i>Enterococcus hirae</i> FDA M19	10541™
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
<i>Bacillus pumilus</i>	14884™
<i>Bordetella bronchiseptica</i>	4617™
<i>Mycobacterium smegmatis</i>	607™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> 3R7089 strain Oxford	9144™
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
<i>Escherichia coli</i> 397E	9637™
<i>Saccharomyces cerevisiae</i>	9763™
EP 5.1.2 BIOLOGICAL INDICATORS OF STERILIZATION	
<i>Bacillus pumilus</i> E601	27142™
<i>Geobacillus stearothermophilus</i>	7953™
<i>Bacillus atrophaeus</i> NRS 1221A	9372™
EP 5.1.3 EFFICACY OF ANTIMICROBIAL PRESERVATION	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
<i>Escherichia coli</i> Crooks	8739™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Candida albicans</i> 3147	10231™
EP MONOGRAPHS: PENICILLAMINE	
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
EP MONOGRAPHS: WATER FOR INJECTIONS	
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
EP MONOGRAPHS: WATER FOR PREPARATION OF EXTRACTS	
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
EP MONOGRAPHS: WATER, HIGHLY PURIFIED	
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
EP MONOGRAPHS: WATER, PURIFIED	
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™

Official Microbial Assays: Japanese Pharmacopeia (JP)

Application	ATCC® No.
JP 16 4.05. I MICROBIOLOGICAL EXAMINATION OF NON-STERILE PRODUCTS: TOTAL VIABLE AEROBIC COUNT	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
JP 16 4.05. II MICROBIOLOGICAL EXAMINATION OF NON-STERILE PRODUCTS: TESTS FOR SPECIFIED MICRO-ORGANISMS	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Escherichia coli</i> Crooks	8739™
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium CDC 6516-60	14028™
<i>Candida albicans</i> 3147	10231™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™ or 19404™
JP 16 4.06 STERILITY TEST	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™ or 19404™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLRI 034	16404™
JP 16 5.02 MICROBIAL LIMIT TEST FOR CRUDE DRUGS	
<i>Escherichia coli</i> Crooks	8739™
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Candida albicans</i> 3147 or 132	10231™ or 2091™
JP 16 MONOGRAPHS: AMPHOTERICIN B	
<i>Saccharomyces cerevisiae</i> NRRL Y-567	9763™
JP 16 MONOGRAPHS: AMPICILLIN SODIUM	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: ARBEKACIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: BACITRACIN	
<i>Micrococcus luteus</i> 130.21	10240™
JP 16 MONOGRAPHS: BEKANAMYCIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: BENZYL PENICILLIN POTASSIUM	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: BLEOMYCIN HYDROCHLORIDE	
<i>Mycobacterium smegmatis</i>	607™
JP 16 MONOGRAPHS: BLEOMYCIN SULFATE	
<i>Mycobacterium smegmatis</i>	607™
JP 16 MONOGRAPHS: CLOXACILLIN SODIUM HYDRATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: CYCLOSERINE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: DIBEKACIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: DICLOXACILLIN SODIUM HYDRATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™

Official Microbial Assays: Japanese Pharmacopeia (JP)

Application	ATCC® No.
JP 16 MONOGRAPHS: ENVIOMYCIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: ERYTHROMYCIN	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: ERYTHROMYCIN ETHYLSUCCINATE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: ERYTHROMYCIN LACTOBIONATE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: ERYTHROMYCIN STEARATE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: FRADIOMYCIN SULFATE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: GENTAMICIN SULFATE	
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™
JP 16 MONOGRAPHS: GRAMICIDIN	
<i>Enterococcus hirae</i> FDA M19	10541™
JP 16 MONOGRAPHS: JOSAMYCIN	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: JOSAMYCIN PROPIONATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: KANAMYCIN MONOSULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: KANAMYCIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: KITASAMYCIN	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: KITASAMYCIN ACETATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: KITASAMYCIN TARTRATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: MICRONOMICIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: MIDECAMYCIN	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: MIDECAMYCIN ACETATE	
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
JP 16 MONOGRAPHS: NYSTATIN	
<i>Saccharomyces cerevisiae</i> NRRL Y-567	9763™
JP 16 MONOGRAPHS: PEPLOMYCIN SULFATE	
<i>Mycobacterium smegmatis</i>	607™
JP 16 MONOGRAPHS: RIBOSTAMYCIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: ROKITAMYCIN	
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
JP 16 MONOGRAPHS: SODIUM FUSIDATE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
JP 16 MONOGRAPHS: SPECTINOMYCIN HYDROCHLORIDE HYDRATE	
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> PCI 602	10031™

Official Microbial Assays: Japanese Pharmacopeia (JP)

Application	ATCC® No.
JP 16 MONOGRAPHS: SPIRAMYCIN ACETATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: STREPTOMYCIN SULFATE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: SULBENICILLIN SODIUM	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: TEICOPLANIN	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: TOBRAMYCIN	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: VANCOMYCIN HYDROCHLORIDE	
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™
JP 16 MONOGRAPHS: ZINOSTATIN STIMALAMER	
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™
JP 16 G3 BIOTECHNOLOGICAL/BIOLOGICAL PRODUCTS: MYCOPLASMA TESTING FOR CELL SUBSTRATES USED FOR THE PRODUCTION OF BIOTECHNOLOGICAL/BIOLOGICAL PRODUCTS	
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	15531™
<i>Mycoplasma orale</i> CH 19299	23714™
<i>Mycoplasma hyorhinis</i> DBS 1050 or BTS-7	29052™ or 17981™
JP 16 G4 MICROORGANISMS: PRESERVATIVES-EFFECTIVENESS TESTS	
<i>Escherichia coli</i> Crooks	8739™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Candida albicans</i> 3147	10231™
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™
JP 16 G4 MICROORGANISMS: TERMINAL STERILIZATION AND STERILIZATION INDICATORS	
<i>Brevundimonas diminuta</i> FDA strain PCI 818	19146™
<i>Geobacillus stearothermophilus</i>	7953™
<i>Bacillus atrophaeus</i> NRS 1221A	9372™
JP 16 G8 QUALITY CONTROL OF WATER FOR PHARMACEUTICAL USE	
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™
<i>Escherichia coli</i> Crooks	8739™
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™
<i>Pseudomonas protegens</i> 1	17386™
<i>Methylobacterium extorquens</i> NBRC 15911	BAA-2500™

STRAINS LISTED BY GENUS AND SPECIES

Strains for Pharmaceutical Microbiology

Strain	ATCC® No.	Application
<i>Acholeplasma laidlawii</i> PG8	23206™	EP 2.6.7 Mycoplasmas, USP Chapter <63>, Mycoplasma test
<i>Aspergillus brasiliensis</i> WLRI 034(120)	16404™	EP 5.1.3 Efficacy of antimicrobial preservation; EP 2.6.1 Sterility; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.12 Microbial examination of non-sterile products: Microbial Enumeration Tests; USP 37 Disinfectants and Antiseptics <1072> ; USP 37 Microbial limit tests – Dietary supplements - Preparatory testing, Growth Promotion Testing <2021>; <51> Antimicrobial Effectiveness Testing, USP 37-NF32, 2014.; USP 37, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests <61>; Sterility Tests <71> USP 37-NF32, 2014; JP 16 4.05 I Microbiological Examination of Non-sterile Products: Total Viable Aerobic Count; JP 16 4.06 Sterility test; JP 16 G4 Microorganisms: Preservatives-Effectiveness Test
<i>Bacillus atrophaeus</i> NRS 1221A	9372™	EP 5.1.2 Biological indicators of sterilization; USP Chapter <1211> Sterilization and sterility assurance of compendial articles, USP 37-NF32, 2014; USP Chapter <1035> Biological Indicators, USP 37-NF32, 2014; USP Monographs: Biological Indicator for Dry-Heat Sterilization, Paper Carrier, USP 37-NF32, 2014; USP Monographs: Biological Indicator for Ethylene Oxide Sterilization, Paper Carrier, USP 37-NF32, 2014; JP 16 G4 Microorganisms: Terminal Sterilization and Sterilization Indicators
<i>Bacillus pumilus</i> E601	27142™	EP 5.1.2 Biological indicators of sterilization; USP Chapter<1211>, Sterilization and sterility assurance of compendial articles
<i>Bacillus pumilus</i>	14884™	EP 2.7.2 Microbiological assays of antibiotics
<i>Bacillus subtilis</i> subsp. <i>spizizenii</i> NRS 231	6633™	EP 2.6.1 Sterility; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.31 Microbiological Examination of Herbal Medicinal Products for Oral Use; EP 2.7.2 Microbiological assays of antibiotics; EP 2.6.12 Microbial examination of non-sterile products: Microbial Enumeration Tests; EP 2.6.13 Microbial examination; EP Monographs: Water for Injections; EP Monographs: Water for Preparation of Extracts; EP Monographs: Water, Highly Purified; EP Monographs: Water, Purified; USP 37 Microbial limit tests – Dietary supplements - Preparatory testing, Growth Promotion Testing <2021>; <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; USP 37, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests <61>; <71> Sterility Tests USP 37-NF32, 2014; JP 16 4.05 I. Microbiological examination of non-sterile products: total viable aerobic count; JP 16 4.06 Sterility test; JP 16 5.02 Microbial limit test for crude drugs; JP 16: Ampicillin sodium, Arbekacin sulfate, Bekanamycin sulfate, Cloxacillin sodium hydrate, Cycloserine, Dibekacin sulfate, Dicloxacillin sodium hydrate, Enviomycin sulfate, Josamycin, Josamycin propionate, Kanamycin monosulfate, Kanamycin sulfate, Kitasamycin, Kitasamycin acetate, Kitasamycin tartrate, Micronomicin sulfate, Midecamycin, Ribostamycin sulfate, Spiramycin acetate, Streptomycin sulfate, Sulbenicillin sodium, Teicoplanin, Tobramycin, Vancomycin hydrochloride
<i>Bacillus subtilis</i> PRD 66	19659™	USP Chapter <1072> Disinfectants and Antiseptics
<i>Bacteroides fragilis</i> VPI 2553	25285™	EP 2.6.27 Microbiological Control of Cellular Products
<i>Bacteroides vulgatus</i>	8482™	USP Chapter <71> Sterility Tests USP 37-NF32, 2014
<i>Bordetella bronchiseptica</i> NRRL B-140	4617™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014
<i>Brevundimonas diminuta</i> FDA strain PCI 818	19146™	USP Chapter <1211> Sterilization and sterility assurance of compendial articles. USP 37-NF32, 2014; JP 16 G4 Microorganisms: Terminal Sterilization and Sterilization Indicators
<i>Candida albicans</i> 132	2091™	EP 2.6.1 Sterility; USP Chapter <1072> Disinfectants and Antiseptics; JP 16 5.02 Microbial Limit Test for Crude Drugs
<i>Candida albicans</i> 3147	10231™	EP 5.1.3 Efficacy of antimicrobial preservation; EP 2.6.1 Sterility; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.12 Microbial examination of non-sterile products: Microbial Enumeration Tests; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; USP Chapter <1072> Disinfectants and Antiseptics ; USP 37 Microbial limit tests – Dietary supplements - Preparatory testing, Growth Promotion Testing <2021>; <51> Antimicrobial Effectiveness Testing, USP 37-NF32, 2014. ; USP 37, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests <61>; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; Sterility Tests <71> USP 37-NF32, 2014; JP 16 4.05 I. Microbiological Examination of Non-Sterile Products: Total Viable Aerobic Count; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-organisms; JP 16 4.06 Sterility Test; JP 16 5.02 Microbial Limit Test for Crude Drugs; JP 16 G4 Microorganisms: Preservatives-Effectiveness Tests
<i>Clostridium sporogenes</i> L.S. McClung 2006	11437™	EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; USP 37, Chondroitin Sulfate Sodium; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; Sterility Tests <71> USP 37-NF32, 2014; Dietary Supplements: Chondroitin Sulfate Sodium, USP 37-NF32, 2014; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-organisms; JP 16 4.06 Sterility Test

Strains for Pharmaceutical Microbiology

Strain	ATCC® No.	Application
<i>Clostridium sporogenes</i>	19404™	EP 2.6.1 Sterility; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; USP 37, Chondroitin Sulfate Sodium; USP 37, Microbiological Examination of Nonsterile Products: Test; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; <71> Sterility Tests USP 37-NF32, 2014; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-Organisms; JP 16 4.06 Sterility Test
<i>Enterococcus hirae</i> FDA M19	10541™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Gramicidin
<i>Escherichia coli</i> MacLeod	10536™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; Tetracycline Hydrochloride and Novobiocin Sodium Tablets. USP 37-NF32, 2014
<i>Escherichia coli</i> AMC 198	11229™	USP Chapter <1072> Disinfectants and Antiseptics USP 37- NF32, 2014
<i>Escherichia coli</i> Crooks	8739™	EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; EP 2.6.31 Microbiological Examination of Herbal Medicinal Products for Oral Use; EP 5.1.3 Efficacy of antimicrobial preservation; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; USP 37 Microbial limit tests – Dietary supplements - Preparatory testing, Growth Promotion Testing <2021>; <51> Antimicrobial Effectiveness Testing, USP 37-NF32, 2014; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-Organisms; JP 16 5.02 Microbial Limit Test for Crude Drugs; JP 16 G4 Microorganisms: Preservatives-Effectiveness Tests; JP 16 G8 Quality Control of Water for Pharmaceutical Use
<i>Escherichia coli</i> 397E	9637™	EP 2.7.2 Microbiological assays of antibiotics
<i>Geobacillus stearothermophilus</i>	7953™	EP 5.1.2 Biological indicators of sterilization; USP Chapter <1211> Sterilization and sterility assurance of compendial articles, USP 37-NF32, 2014; USP Monographs: Biological Indicators for Steam Sterilization, Paper Carrier, USP 37-NF32, 2014; USP Chapter <1229.1> Steam Sterilization by Direct Contact, USP 37-NF32, 2014; JP 16 G4 Microorganisms: Terminal Sterilization and Sterilization Indicators
<i>Geobacillus stearothermophilus</i> NCA 26	12980™	USP Chapter <1229.1> Steam Sterilization by Direct Contact, USP 37-NF32, 2014; USP Chapter <1035> Biological Indicators, USP 37-NF32, 2014; USP Monographs Biological Indicator for Steam Sterilization, Paper carrier, USP37-NF32, 2014
Hybrid Moloney/Amphotropic murine leukemia virus (Mo/A-MuLv) 4070A envelope strain	VR-1448™	USP Chapter:<1046> Cell and Gene Therapy Products - Analytical Methodologies, USP 37-NF32, 2014
Hybrid Moloney/Amphotropic murine leukemia virus (Mo/A-MuLv) 4070A envelope strain	VR-1450™	USP Chapter:<1046> Cell and Gene Therapy Products - Analytical Methodologies, USP 37-NF32, 2014
<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> PCI 602	10031™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Spectinomycin Hydrochloride Hydrate
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; Penicillamine USP 37-NF32, 2014; <71> Sterility Tests USP 37-NF32, 2014
<i>Lactobacillus leichmannii</i> 313	7830™	USP Chapter <171> Vitamin B12 Activity Assay, USP 37-NF32, 2014; USP Dietary Supplements: Oil-and Water-Soluble Vitamins with Minerals Tablets, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Tablets, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins Tablets, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Tablets, USP 37-NF32, 2014
<i>Lactobacillus plantarum</i> 17-5	8014™	USP Chapter <91> Calcium Pantothenate Assay, USP 37-NF32, 2014; USP Chapter <441> Niacin or Niacinamide Assay. USP 37-NF32, 2014; Dietary Supplements: Oil-and Water-Soluble Vitamins with Minerals Tablets, USP 37-NF32, 2014.; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Tablets, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins Tablets, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Tablets, USP 37-NF32, 2014
<i>Methylobacterium extorquens</i> NBRC 15911	BAA-2500™	JP 16 G8 Quality Control of Water for Pharmaceutical Use
<i>Micrococcus luteus</i> 130.21	10240™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Bacitracin
<i>Kocuria rhizophila</i> FDA strain PCI 1001	9341™	JP 16: Midecamycin Acetate, Rokitamycin, Zinostatin Stimalamer

Strains for Pharmaceutical Microbiology

Strain	ATCC® No.	Application
<i>Mycobacterium smegmatis</i>	607™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Bleomycin Hydrochloride, Bleomycin Sulfate, Peplomycin Sulfate
<i>Mycoplasma fermentans</i> PG 18	19989™	USP Chapter <63> Mycoplasma test, EP 2.6.7 Mycoplasmas
<i>Mycoplasma gallisepticum</i> PG31	19610™	EP 2.6.7 Mycoplasmas; USP Chapter <63> Mycoplasma test
<i>Mycoplasma hyorhinis</i> BTS-7	17981™	EP 2.6.7 Mycoplasmas; USP Chapter <63> Mycoplasma test; JP 16 G3 Biotechnological/Biological Products: Mycoplasma Testing for Cell Substrates used for the Production of Biotechnological/Biological Products
<i>Mycoplasma hyorhinis</i> DBS 1050	29052™	JP 16 G3 Biotechnological/Biological Products: Mycoplasma Testing for Cell Substrates used for the Production of Biotechnological/Biological Products
<i>Mycoplasma orale</i> CH 19299	23714™	EP 2.6.7 Mycoplasmas; USP Chapter <63> Mycoplasma test; JP 16 G3 Biotechnological/Biological Products: Mycoplasma Testing for Cell Substrates used for the Production of Biotechnological/Biological Products
<i>Mycoplasma pneumoniae</i> FH strain of Eaton Agent	15531™	EP 2.6.7 Mycoplasmas; USP Chapter <63> Mycoplasma test; JP 16 G3 Biotechnological/Biological Products: Mycoplasma Testing for Cell Substrates used for the Production of Biotechnological/Biological Products
<i>Mycoplasma synoviae</i> WVU 1853	25204™	EP 2.6.7 Mycoplasmas; USP Chapter <63> Mycoplasma test
<i>Pediococcus acidilactici</i>	8042™	USP Chapter <115> Dexpanthenol Assay, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Oil- and Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Capsules, USP 37-NF32, 2014; USP Dietary Supplements: Water-Soluble Vitamins with Minerals Oral Solution, USP 37-NF32, 2014; Dietary Supplements: Water-soluble Vitamins Capsules, USP 37-NF32, 2014
<i>Penicillium chrysogenum</i> Wis. 49-133	11709™	USP Chapter <1072> Disinfectants and Antiseptics
<i>Propionibacterium acnes</i> Gerath	11827™	EP 2.6.27 Microbiological Control of Cellular Products
<i>Pseudomonas aeruginosa</i> PRD-10	15442™	USP Chapter <1072> Disinfectants and Antiseptics
<i>Pseudomonas aeruginosa</i>	25619™	USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014
<i>Pseudomonas aeruginosa</i> R. Hugh 813	9027™	EP 5.1.3 Efficacy of antimicrobial preservation; EP 2.6.1 Sterility; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; EP 2.6.31 Microbiological Examination of Herbal Medicinal Products for Oral Use; EP Monographs: Water for Injections; EP Monographs: Water for Preparation of Extracts; EP Monographs: Water, Highly Purified; EP Monographs: Water, Purified; USP 37, Chondroitin Sulfate Sodium; USP, Microbiological Examination of Nonsterile Products: Test; EP 2.6.12 Microbial examination of non-sterile products: Microbial Enumeration Tests; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; USP, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; USP 37 Microbial limit tests – Dietary supplements – Preparatory testing, Growth Promotion Testing <2021>; <51> Antimicrobial Effectiveness Testing, USP 37-NF32, 2014; USP 37, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests <61>; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; <71> Sterility Tests USP 37-NF32, 2014; JP 16 4.05 I. Microbiological Examination of Non-Sterile Products: Total Viable Aerobic Count; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-Organisms; JP 16 4.06 Sterility Test; JP 16 G4 Microorganisms: Preservatives-Effectiveness Tests; JP 16 G8 Quality Control of Water for Pharmaceutical Use
<i>Pseudomonas protegens</i> 1	17386™	JP 16 G8 Quality Control of Water for Pharmaceutical Use
<i>Saccharomyces cerevisiae</i>	2601™	USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014
<i>Saccharomyces cerevisiae</i> NRRL Y-567	9763™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Amphotericin B, Nystatin
<i>Salmonella enterica</i> subsp. <i>enterica</i>	13311™	USP Microbial limit tests – Dietary supplements – Preparatory testing, Growth Promotion Testing <2021>
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Typhimurium CDC 6516-60	14028™	EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; EP 2.6.31 Microbiological Examination of Herbal Medicinal Products; USP Chapter <62>, Microbiological Examination of Nonsterile Products; JP16 4.05 II. Microbiological Examination of Nonsterile Products.
<i>Salmonella enterica</i> subsp. <i>enterica</i> serovar Abony	BAA-2162™	USP Chapter <62> Microbiological examination of non-sterile products: tests for specified microorganisms
<i>Serratia marcescens</i> PCI 1107	14756™	USP Chapter <1211> Sterilization and sterility assurance of compendial articles. USP 37-NF32, 2014
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> PCI 1209/N	12692™	Penicillin G Procaine and Novobiocin Sodium Intramammary Infusion. USP 37-NF32, 2014
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA	29737™	USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; USP 37, Tetracycline Hydrochloride and Novobiocin Sodium Tablets
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> 3R7089 strain Oxford	9144™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014

Strains for Pharmaceutical Microbiology

Strain	ATCC® No.	Application
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> FDA 209	6538™	EP 2.6.1 Sterility testing; EP 2.6.27 Microbiological Control of Cellular Products; EP 2.7.2 Microbiological assays of antibiotics; EP 2.6.12 Microbial examination of non-sterile products: Microbial Enumeration Tests; EP 2.6.13 Microbial examination of non-sterile products: Test for Specified Micro-organisms; EP 2.6.31 Microbiological Examination of Herbal Medicinal Products for Oral Use; USP 37 Disinfectants and Antiseptics <1072> ; USP 37 Microbial limit tests – Dietary supplements - Preparatory testing, Growth Promotion Testing <2021>; <51> Antimicrobial Effectiveness Testing, USP 37-NF32, 2014.; USP 37, Microbiological Examination of Nonsterile Products: Microbial Enumeration Tests <61>; USP 37, Microbiological Examination of Nonsterile Products: Tests for Specified Microorganisms <62>; <71> Sterility Tests USP 37-NF32, 2014; JP 16 4.05 I. Microbiological Examination of Non-Sterile Products: Total Viable Aerobic Count; JP 16 4.05 II. Microbiological Examination of Non-Sterile Products: Tests for Specified Micro-Organisms; JP 16 4.06 Sterility Test; JP 16 5.02 Microbial Limit Test for Crude Drugs, JP 16: Benzylpenicillin Potassium, Erythromycin, Erythromycin Ethylsuccinate, Erythromycin Lactobionate, Erythromycin Stearate, Fradiomycin Sulfate, Sodium fusidate; JP 16 G4 Microorganisms: Preservatives-Effectiveness Tests; JP 16 G8 Quality Control of Water for Pharmaceutical Use
<i>Staphylococcus epidermidis</i> FDA strain PCI 1200	12228™	EP 2.7.2 Microbiological assays of antibiotics; USP Chapter <81> Antibiotics – Microbial Assays, USP 37-NF32, 2014; JP 16: Gentamycin Sulfate
<i>Streptococcus pyogenes</i> Bruno	19615™	EP 2.6.27 Microbiological Control of Cellular Products
<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i> 33114	9610™	EP 2.6.27 Microbiological Control of Cellular Products

Please visit our website www.atcc.org/pharmaceutical for more resources to help you in your pharmaceutical-related research applications.

ATCC QUALITY CONTROL STANDARDS

Search our quality control page at www.atcc.org/qc for strains that are organized by procedural standards and commercial identification systems.

- Standards from organizations and agencies such as CLSI, ASTM International, USP, USDA, FDA, AOAC, and BSI-Global
- Commercial identification systems from companies including bioMérieux, Biolog, Siemens, DuPont Nutrition and Health, and others

ATCC RELATED PRODUCTS

Value-added microbial reagents, such as our microbial freeze media and supplements, as well as nucleic acids offer convenience while saving time and money.

- Mycoplasma Growth Supplement (ATCC® No. 20-2206 and 20-2207) and TSB with 10% glycerol (ATCC® No. 20-2200)
- Genomic DNA and RNA from well-characterized and authenticated microbial strains



Convenient Tools for Quality Control

Developed by the leaders in microbial cultivation and preservation, ATCC® Minis provide a convenient, ready-to-use solution for handling quality control strains. Each glass-free tube contains a ready-to-use strain in glycerol stock, with a 2D barcode for easy tracking and peel-off labels for reliable recordkeeping. It is easy to ensure the quality of your products with ATCC® Minis – just open, plate, and go! Visit us online at www.atcc.org/minis to browse our collection.



Find additional resources at www.atcc.org/pharmaceutical.

BA-0317-07

© 2017 American Type Culture Collection. The ATCC trademark and trade name, any and all ATCC catalog numbers and any other trademarks listed in this publication are trademarks of the American Type Culture Collection unless indicated otherwise.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services or to manufacture commercial products without prior ATCC written approval.

Tel 800.638.6597
703.365.2700
Fax 703.365.2750
Email sales@atcc.org
Web www.atcc.org

Or contact your local distributor

