

Antibiotic Susceptibility of *Pseudomonas aeruginosa* ATCC® BAA-2113™

Antimicrobial	MIC ^a	Interpretation ^b	Antimicrobial	MIC	Interpretation
Amikacin	16	S	Imipenem	0.5	S
Aztreonam	>16	I	Levofloxacin	1	S
Cefepime	4	S	Meropenem	1	S
Ceftazidime	4	S	Piperacillin	16	S
Ciprofloxacin	≤0.5	S	Piperacillin-tazobactam	8/4	S
Doripenem	≤0.5	S	Tobramycin	≤2	I

^a Antibiotic susceptibility was obtained using the Sensititre® brand of antimicrobial titration analysis system.

^b Parameter set: MIC Interpretation Guideline: CLSI M100 ED33:2023