

**ATCC medium: 1432 Potato P-YE thermus medium**

Potato Extract (see below).....200.0 ml  
Peptone.....5.0 g  
Yeast extract.....0.2 g  
Agar (if needed).....20.0 g  
Distilled water.....800.0 ml

Adjust medium for final pH 7.8. Autoclave at 121C for 15 minutes.

*Potato Extract:*

To prepare potato extract, boil 300 g peeled, diced potatoes in 500 ml distilled water for half an hour. Decant the infusion and use in the above medium.