

ATCC medium: 1794 PGA medium

Peptone.....5.0 g
Glucose.....5.0 g
Na₂HPO₄1.2 g
KH₂PO₄0.5 g
NaCl.....3.0 g
Agar.....18.0 g
Potato Decoction (see below).....1.0 L

Adjust medium for final pH 7.0. Autoclave at 121C for 15 minutes.

Potato Decoction:

Boil 200 g peeled, diced potatoes in 1.0 L distilled water for 20 minutes. Filter through cheesecloth and return volume of filtrate to 1.0 L.