

ATCC medium: 1890 LB medium (ATCC medium 1065) with 12 mcg/ml tetracycline and 100 mcg/ml ampicillin

ATCC Medium 1065 (see below) with 12.0 mcg/ml tetracycline and 100.0 mcg/ml ampicillin

Prepare a filter-sterilized solutions of tetracycline and ampicillin and aseptically; add enough of each to the sterile basal medium to achieve the stated concentrations.

ATCC Medium 1065:

Tryptone (BD 211705).....	10.0 g
Yeast Extract (BD 212750).....	5.0 g
NaCl.....	10.0 g
Distilled water.....	1.0 L

Autoclave at 121C for 15 minutes.