

ATCC medium: 2634 Haloterrigena medium

Yeast extract.....	5.0 g
Casamino acids.....	5.0 g
Na-glutamate.....	1.0 g
KCl.....	2.0 g
Na ₃ -citrate.....	3.0 g
MgSO ₄ . 7H ₂ O.....	20.0 g
NaCl.....	200.0 g
FeCl ₂ . 4H ₂ O.....	36.0 g
MnCl ₂ . 4H ₂ O.....	0.4 g
Agar.....	20.0 g
Distilled water.....	1.0 L

Adjust the pH to 7.0 to 7.2. Autoclave at 121°C for 15 minutes.