

**ATCC Medium: 2761 Natronolimnibus Medium**

KH <sub>2</sub> PO <sub>4</sub> .....	0.3 g
KCl.....	1.0 g
NaCl.....	200 g
Yeast Extract.....	0.5 g
Tryptone.....	0.5 g
Casamino Acids.....	0.5 g
Na <sub>2</sub> CO <sub>3</sub> .....	6.8 g
NaHCO <sub>3</sub> .....	3.8 g
DI Water.....	1000 ml

Autoclave medium at neutral pH at 121°C.

Aseptically add the following stock solutions:

Sodium Pyruvate (20% wt/vol.).....	2 ml
Glucose (20% wt/vol.).....	4 ml

Adjust pH to 9.5-10.5 after medium has cooled.