

## **ATCC Medium 2901: Nutrient Medium + 25 $\mu$ g/mL Imipenem**

### **Agar Medium**

Nutrient Agar (BD 213000) .....	23.0 g
DI Water.....	950.0 mL
Imipenem Solution (See below).....	50.0 mL

Autoclave at 121° C. Cool to ~50°C and aseptically add Imipenem Solution.  
Dispense as required.

### **Broth Medium**

Nutrient Broth (BD 234000) .....	8.0 g
DI Water.....	950.0 mL
Imipenem Solution (See below).....	50.0 mL

Autoclave at 121° C. Cool to ~50°C and aseptically add Imipenem Solution.  
Dispense as required.

### **Imipenem Solution (25 $\mu$ g/mL)**

Imipenem.....	25.0 mg
DI Water .....	50.0 mL

Filter sterilize.

\*This product is commercially available from BD.