

**ATCC medium: 564 GYPT medium**

Glucose.....5.0 g  
Yeast extract.....6.0 g  
Peptone, Bacto (BD 211677)...2.0 g  
HEPES (Sigma H 3375).....4.0 g  
Sodium thioglycollate.....0.5 g  
Agar.....8.0 g  
Distilled water.....1.0 L

Adjust pH to 7.4 with KOH. Tube under 97% N<sub>2</sub>, 3% H<sub>2</sub>, cap with butyl rubber stoppers. Autoclave in press for 15 minutes under fast exhaust.