

ATCC medium: 653 Modified trypticase soy (ATCC medium 646) with 10% horse serum

Aseptically add 100 ml sterile horse serum to ATCC Medium 646 (see below)

ATCC Medium 646:

Trypticase Peptone (BD 211921)	17.0 g
Phytone Peptone (BD 211906)	3.0 g
NaCl	5.0 g
K ₂ HPO ₄	2.5 g
Yeast extract	4.0 g
Agar	15.0 g
Distilled water	1.0 L

Autoclave at 121C for 15 minutes.