

**ATCC medium: 89 *Micrococcus-Sarcina* medium**

K <sub>2</sub> HPO <sub>4</sub> .....	0.2 g
Sodium citrate . 2H <sub>2</sub> O .....	0.5 g
Sodium acetate . 3H <sub>2</sub> O .....	0.3 g
Sodium succinate . 6H <sub>2</sub> O .....	2.0 g
Trypticase Peptone (BD 211921) .....	5.0 g
Yeast autolysate .....	1.0 g
Starch (laundry) .....	2.0 g
Agar .....	16.0 g
Distilled water .....	1.0 L

Adjust pH to 7.0.  
Autoclave at 121C for 15 minutes.