

ATCC Medium: 2553 Tryptic Soy w/ 150mcg/ml of Ampicillin and 150mcg/ml Streptomycin

For Broth Tryptic Soy Broth (Difco).....30.0g
 Distilled water.....900 mls
 Mix thoroughly and autoclave @ 121 C

For Agar: Tryptic Soy Agar (Difco).....40.0g
 Distilled water.....900 mls
 Heat with frequent agitation and boil for 1 minute. Autoclave @ 121°C.

Aseptically add to sterile media:

Ampicillin.....0.15g
Streptomycin.....0.15g
Distilled water.....100 mls

Filter Sterilize

Dispense complete media into appropriate vessel.