

ATCC Medium 5129: Haemophilus Test Medium (broth)

Solution A: Bovine Hematin (*dissolved in NaOH*), 15.0 mg
Mueller Hinton Broth (Difco 275730), 21.0 g
Distilled Water, 900 ml

Sterilize at 121°C for 15 minutes. Cool to 50°C.

Solution B: NAD, 15.0 mg
Yeast Extract, 5.0 g
Distilled Water, 100 ml

Filter sterilize solution B and add to cooled solution A.

Adjust pH for final 7.3 +/- 0.1. Dispense.