
## Technical Data Sheet: SUP-T1 NFAT-Luc2


ATCC® Number	CRL-1942-NFAT-LUC2™
Organism	Homo sapiens
Tissue/Disease Source	T-cell; Lymphoblastic Lymphoma
Product Description	The SUP-T1 cell line (ATCC CRL-1942) is commonly used in immuno-oncology research and endogenously expresses a high level of programmed cell death protein 1 (PD-1). This luciferase reporter cell line was derived from the parental line CRL-1942 by stably expressing the firefly luciferase gene (luc2) under control of a nuclear factor of activated T-cells (NFAT) promoter. This cell line was established through lentiviral transduction and single cell cloning. The cells, upon stimulation with PMA and ionomycin, express high levels of enzymatically active luciferase protein, which can be detected via in vitro bioluminescence assays. This reporter cell line is useful for monitoring the activity of NFAT signaling pathways that regulate a wide range of cell responses including immune cell activation.
Application	Enabling sensitive and quantitative assessment of signal transduction makes this reporter cell line ideal for in vitro bioluminescence assays to study immune response in cell lines overexpressing PD-1, development of new drugs, and safety evaluation of new chemicals and drugs.

## In vitro activation of luciferase expression by PMA and lonomycin




**Figure 1. In vitro activation of luciferase expression by PMA and lonomycin.** Luciferase expression from SUP-T1 cells upon signaling activation by stimulation with 50 ng/mL of PMA and 10 μg/mL of lonomycin to activate luciferase expression after 6 hours incubation to demonstrate response sensitivity in NFAT response element. N=3 in all experiments.

## **Expression of PD-1**



**Figure 2. Expression of PD-1.** Flow cytometry analysis was performed to assess the antigen expression levels of PD-1 (pink) on the tumor cell lines compared to isotype controls (blue).

## **Cell morphology**



**Figure 3: Cell morphology of SUP-T1 parental and SUP-T1 NFAT-Luc2.** Cells were maintained in ATCC recommended culture conditions. Cell morphology was observed under microscopy and images were captured by digital camera.

© 2024 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are owned by the American Type Culture Collection unless indicated otherwise.

These products are for laboratory use only. Not for human or diagnostic use. ATCC products may not be resold, modified for resale, used to provide commercial services, or to manufacture commercial products without prior ATCC written approval.