

ATCC Medium: 2836 BHI Agar/Broth w/ 100 mcg/ml Rifampicin

Agar Medium

Brain Heart Infusion Agar (BD 211065).....	52 g
Rifampicin Solution (see below, add aseptically)....	50 ml
DI Water.....	950 ml

Autoclave at 121°C. Cool sterilized medium to ~50°C. Aseptically add 50 ml of filter sterilized rifampicin solution. Gently mix and dispense as required.

Broth Medium

Brain Heart Infusion (BD 237500).....	37.0 g
Rifampicin Solution (see below, add aseptically)....	50 ml
DI Water.....	950 ml

Autoclave at 121°C. Cool sterilized medium to ~50°C. Aseptically add 50 ml of filter sterilized rifampicin solution. Gently mix and dispense as required.

Rifampicin Solution

Rifampicin.....	0.1 g
DI Water.....	50 ml

Filter sterilize and add to sterile base.