

ATCC medium: 1543 Mineral medium with santonin

K ₂ HPO ₄	6.3 g
KH ₂ PO ₄	1.82 g
NH ₄ NO ₃	1.0 g
MgSO ₄ . 7H ₂ O	0.1 g
CaCl ₂ . 2H ₂ O	0.75 g
FeSO ₄ . 7H ₂ O	0.06 g
MnSO ₄ (anhydrous)	600.0 mcg
Na ₂ MoO ₄ . 2H ₂ O.....	600.0 mcg
alpha-santonin.....	4.0 g
Distilled water.....	1.0 L

Adjust to pH 7.0. Autoclave at 121C for 15 minutes.