

ATCC medium: 1586 1-Chlorobutane medium

KH₂PO₄1.5 g
Na₂HPO₄ . 12H₂O.....1.5 g
NH₄NO₃4.0 g
MgSO₄ . 7H₂O10.0 mg
CaSO₄ . 2H₂O10.0 mg
FeSO₄ . 7H₂O5.0 mg
Yeast extract.....5.0 mg
1-Chlorobutane*.....0.2 %
Distilled water.....1.0 L

Adjust medium for final pH 7.0. Autoclave at 121C for 15 minutes.

*Filter-sterilize separately and add aseptically to sterile basal medium.