

ATCC medium: 1806 CCY Modified medium

Yeast extract.....30.0 g
Casamino acids.....20.0 g
Na₂HPO₄2.48 g
KH₂PO₄0.41 g
MgSO₄ . 7H₂O20.0 mg
MnSO₄ . H₂O7.5 mg
FeSO₄ . 7H₂O6.4 mg
Citric acid.....6.4 mg
Distilled water.....900.0 ml

Adjust for final pH 7.3. Autoclave at 121C for 15 minutes. Cool to room temperature and aseptically add the following filter-sterilized solution:

Sodium pyruvate.....23.2 g
Distilled water.....100.0 ml