

ATCC medium: 1907 Rhizomonas medium plus tetracycline

ATCC 1700 Rhizomonas medium

MgSO₄ . 7H₂O0.5 g
KNO₃0.5 g
K₂HPO₄1.0 g
Glucose.....2.5 g
Trypticase Peptone (BD 211921).....5.0 g
Ca(NO₃)₂ . 4H₂O.....60.0 mg
Agar, Noble (BD 214230).....11.0 g
Distilled water.....1.0 L

Adjust medium to pH of 7.2. Autoclave at 121C for 15 minutes.

To above solution add
Tetracycline.....10.0 mcg

Filter-sterilize.