

**ATCC medium: 1930 Photosynthetic *Rhizobium* Medium**

K<sub>2</sub>HPO<sub>4</sub> .....0.5 g  
MgSO<sub>4</sub> . 7H<sub>2</sub>O .....0.2 g  
NaCl.....0.1 g  
Glutamic Acid.....9.55 g  
Yeast Extract.....1.0 g  
FeSO<sub>4</sub> . 7H<sub>2</sub>O .....0.009 g  
Agar (for plates).....15.0 g  
Agar (for slants).....20.0 g  
Distilled/deionized water.....1.0 L

Adjust to pH 6.9. Autoclave at 121C for 15 minutes.