

ATCC medium: 543 *Rhodopseudomonas* medium

Sodium succinate.....2.5 g
K₂HPO₄0.9 g
KH₂PO₄0.6 g
MgSO₄ . 7H₂O0.2 g
(NH₄)₂SO₄.....1.25 g
CaCl₂0.07 g
Ferric citrate.....3.0 mg
EDTA.....2.0 mg
Yeast extract.....0.5 g
Distilled water.....1.0 L

Adjust pH to 7.0.
Autoclave at 121C for 15 minutes.