

ATCC medium: 80 Lactic bacteria broth

Trypticase Peptone (BD 211921).....10.0 g
Yeast autolysate.....5.0 g
Sodium acetate.....12.0 g
Glucose.....10.0 g
Solution A (see below).....5.0 ml
Solution B (see below).....5.0 ml
Distilled water to.....1.0 L

Adjust pH to 5.1-5.3.
Autoclave at 121C for 15 minutes.

Solution A:

K_2HPO_4 10.0 %
 KH_2PO_4 10.0 %

Solution B:

$MgSO_4$ 4.0 %
 $NaCO_3$ 0.2 %
 $FeSO_4$ 0.2 %
 $MnSO_4$ 0.2 %