

**ATCC Medium: 2714 Middlebrook 7H9 Broth w/ ADC Enrichment and 40mM Sodium Pyruvate**

Middlebrook 7H9 Broth.....	4.7 g
Glycerol.....	2.0 ml
Sodium pyruvate.....	4.4 g
DI Water.....	900.0 ml

Add 4.7 grams of dehydrated Middlebrook 7H9 medium and 4.4 grams Sodium pyruvate to 900 ml of water containing 2.0 ml of glycerol. Mix well and autoclave at 121°C. Cool to 50-55°C.

Aseptically add 100 ml of the following filter sterilized ADC Enrichment:

**ADC Enrichment Formulation**

Bovine Serum Albumin, Frac V.....	5.0 g
Dextrose.....	2.0 g
Catalase(beef).....	3 mg
DI Water.....	100 ml