

## **ATCC Medium: 3 Nutrient Agar/Broth**

### **Agar Medium**

Nutrient Agar (BD 213000).....23 g  
DI Water.....1000 ml

Autoclave at 121°C.

### **Broth Medium**

Nutrient Broth (BD cat 234000).....8.0 g  
DI Water.....1000 ml

Autoclave at 121°C.

\*This product is commercially available from BD.

To make medium from scratch, follow formulation below:

### **Nutrient Agar Composition**

Beef Extract.....3.0 g  
Peptone.....5.0 g  
Agar.....15.0 g

Final pH 6.8 +/- 0.2.

\*Omit agar for broth medium.