

ATCC Medium: 2854 Nutrient Agar/Broth with 300 mcg/ml Ampicillin

Agar Medium

Nutrient Agar (BD 213000).....23 g
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

Ampicillin Solution

Ampicillin.....300 mg
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**

Broth Medium

Nutrient Broth (BD 234000).....8.0 g
DI Water.....950 ml

Autoclave at 121°C and allow to cool down to 45-50°C. Aseptically add the following filter sterilized solution:

Ampicillin Solution

Ampicillin.....300 mg
DI Water.....50 ml

Dispense as required. **Store at 4°C. Do not store longer than 7-10 days.**